



09D0.HC1 Hand Held Crimping Tool

A compact, hand-held design that is easy to use.
Will crimp pumps from 13mm to 20mm. Supplied with 2 collet sizes
Execution of bottle crimping is fast and simple making the tool ideal for small-scale production or lab use.
This crimping tool allows crimping with the perfume bottle pump actuator (push-button) in place, (for most common sizes) simplifying the crimping process since there is no need to remove and replace it.
A simple adjustment allowing the crimp to be used with multiple bottles and pump sizes.
Dimensions: 190mm long, 40mm average diameter, 1.1kg weight.
This tool has helped many small perfume packing businesses to start manufacturing perfume without the expense of a large crimping machine.



09D.0CT Collar Fitting Tool

Low cost tool for fitting collars and crimpless pumps
simple and easy to use.
Just place the bottom part of the tool over the collar or crimpless pump, hold the top on the knob, and use the silver sliding weight to tap the collar down into position.
Weight- 430 g